

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

APPLICATION FOR LETTERS PATENT

Inventors: Marc De Beuckeleer

Title: **METHODS AND KITS FOR IDENTIFYING ELITE EVENT
GAT-ZM1 IN BIOLOGICAL SAMPLES**

EXPRESS MAIL

Mailing Label Number EF308522607US

Date of Deposit January 11, 2001

I hereby certify that this paper or fee is being
deposited with the United States Postal Service

"Express Mail Post Office to Addressee" Service

under 37 CFR 1.10 on the date indicated above and

is addressed to the Commissioner of Patents and

Trademarks, Washington, D.C. 20231

Kenneth Patterson

(Typed or printed name of person

mailing paper or fee)

Kenneth Patterson

(Signature of person mailing paper or fee)

Susan Lehnhardt
Reg. No. 33,943
Frommer Lawrence & Haug, LLP
745 Fifth Avenue
New York, New York 10151
(212) 588-0800

TITLE OF THE INVENTION

Methods and kits for identifying elite event GAT-ZM1 in biological samples.

CROSS-REFERENCE TO RELATED APPLICATIONS

This application is a continuation of U.S. Serial No. 09/481,049 incorporated

5 herein by reference.

STATEMENT REGARDING FEDERALLY SPONSORED RESEARCH AND DEVELOPMENT

Not applicable.

BACKGROUND OF THE INVENTION

10 The phenotypic expression of a transgene in a plant is determined both by the
structure of the gene itself and by its location in the plant genome. At the same time the
presence of the transgene (in a foreign DNA) at different locations in the genome will
influence the overall phenotype of the plant in different ways. The agronomically or
industrially successful introduction of a commercially interesting trait in a plant by genetic
15 manipulation can be a lengthy procedure dependent on different factors. The actual
transformation and regeneration of genetically transformed plants are only the first in a
series of selection steps, which include extensive genetic characterization, breeding, and
evaluation in field trials, eventually leading to the selection of an elite event.

20 The unequivocal identification of an elite event is becoming increasingly important
in view of discussions on Novel Food/Feed, segregation of GMO and non-GMO products
and the identification of proprietary material. Ideally, such identification method is both
quick and simple, without the need for an extensive laboratory set-up. Furthermore, the
method should provide results that allow unequivocal determination of the elite event
without expert interpretation, but which hold up under expert scrutiny if necessary.

GAT-ZM1 was selected as an elite event in the development of corn resistant to the herbicide Liberty®, by transformation of corn with plasmid pUC/Ac comprising the pat gene encoding tolerance to phosphinothricin. It is commercially sold as Liberty Link® maize, such as, for instance, Liberty Link® A6460LL sold by AgriGold/Akin Seed Company. The tools for use in simple and unequivocal methods for identification elite event GAT-ZM1 in biological samples are described herein.

SUMMARY OF THE INVENTION

The present invention relates to methods for identifying elite event GAT-ZM1 in biological samples, which methods are based on primers or probes that specifically recognize the 5' and/or 3' flanking sequence of GAT-ZM1.

More specifically, the invention relates to a method comprising amplifying a sequence of a nucleic acid present in biological samples, using a polymerase chain reaction with at least two primers, one of which recognizes the 5' or 3' flanking region of GAT-ZM1, the other which recognizes a sequence within the foreign DNA, to obtain a DNA fragment of between 100 and 350 bp. Preferably, the primers recognize a sequence within the 5' flanking region of GAT-ZM1, most preferably within the 5' flanking region of SEQ ID NO: 6, and a sequence within the foreign DNA, respectively. Especially preferably, the primer recognizing the 5' flanking region comprises the nucleotide sequence of SEQ ID NO: 11 and the primer recognizing a sequence within the foreign DNA comprises the nucleotide sequence of SEQ ID NO: 12 described herein.

The present invention more specifically relates to a method for identifying elite event GAT-ZM1 in biological samples, which method comprises amplifying a sequence of a nucleic acid present in a biological sample, using a polymerase chain reaction with two

primers having the nucleotide sequence of SEQ ID NO: 11 and SEQ ID NO: 12 respectively, to obtain a DNA fragment of between 180 and 220 bp, preferably of about 200 bp.

The present invention further relates to the specific flanking sequences of GAT-ZM1 described herein, which can be used to develop specific identification methods for GAT-ZM1 in biological samples. More particularly, the invention relates to the 5' and or 3' flanking regions of GAT-ZM1 that can be used for the development of specific primers and probes. The invention further relates to identification methods for the presence of GAT-ZM1 in biological samples based on the use of such specific primers or probes.

The invention further relates to kits for identifying elite event GAT-ZM1 in biological samples, said kits comprising at least one primer or probe which specifically recognizes the 5' or 3' flanking region of GAT-ZM1.

Preferably the kit of the invention comprises, in addition to a primer that specifically recognizes the 5' or 3' flanking region of GAT-ZM1, a second primer that specifically recognizes a sequence within the foreign DNA of GAT-ZM1, for use in a PCR identification protocol. Preferably, the kit of the invention comprises two specific primers, one of which recognizes a sequence within the 5' flanking region of GAT-ZM1, most preferably within 5' flanking region of SEQ ID NO: 6, and the other which recognizes a sequence within the foreign DNA. Especially preferably, the primer recognizing the 5' flanking region comprises the nucleotide sequence of SEQ ID NO: 11 and the primer recognizing the transgene comprises the nucleotide sequence of SEQ ID NO: 12 herein.

The invention further relates to a kit for identifying elite event GAT-ZM1 in biological samples, said kit comprising the PCR primers having the nucleotide sequence of

SEQ ID NO: 11 and SEQ ID NO: 12 for use in the GAT-ZM1 PCR identification protocol described herein.

The invention also relates to a kit for identifying elite event GAT-ZM1 in biological samples, which kit comprises a specific probe having a sequence which
5 corresponds (or is complementary to) a sequence having between 80% and 100% sequence identity with a specific region of GAT-ZM1. Preferably the sequence of the probe corresponds to a specific region comprising part of the 5' or 3' flanking region of GAT-ZM1. Most preferably the specific probe has (or is complementary to) a sequence having between 80% and 100% sequence identity to the sequence between nucleotide 286
10 and 466 of SEQ ID NO: 6.

The methods and kits encompassed by the present invention can be used for different purposes such as, but not limited to, the following: to identify GAT-ZM1 in plants, plant material or in products such as, but not limited to food or feed products (fresh or processed) comprising or derived from plant material; additionally or alternatively, the
15 methods and kits of the present invention can be used to identify transgenic plant material for purposes of segregation between transgenic and non-transgenic material; additionally or alternatively, the methods and kits of the present invention can be used to determine the quality (i.e. percentage pure material) of plant material comprising GAT-ZM1.

The invention further relates to the 5' and/or 3' flanking regions of GAT-ZM1 as
20 well as to the specific primers and probes developed from the 5' and/or 3' flanking sequences of GAT-ZM1.

BRIEF DESCRIPTION OF THE DRAWING

The following Examples, not intended to limit the invention to specific embodiments described, may be understood in conjunction with the accompanying Figure, incorporated herein by reference, in which:

Figure 1 depicts the scoring of unknowns using the PCR Identification protocol developed for GAT-ZM1. Loading sequence of the gel: Unknowns: lanes 1, 2, 5, 6, 8, 11, 13, 14, DNA samples from corn plants comprising the transgenic event GAT-ZM1; lanes 4, 9, 10, 12, DNA samples from a corn plants not comprising elite event GAT-ZM1; lane 3, PCR failure. Control lanes: lanes 19, 21, control DNA samples from corn plants comprising elite event GAT-ZM1; lanes 20, 22, control DNA samples from wild-type corn plants; lane 23, no template control; lane 24, molecular weight marker.

DETAILED DESCRIPTION OF THE INVENTION

The incorporation of a recombinant DNA molecule in the plant genome typically results from transformation of a cell or tissue (or from another genetic manipulation). The particular site of incorporation is either due to “random” integration or is at a predetermined location (if a process of targeted integration is used).

The DNA introduced into the plant genome as a result of transformation of a plant cell or tissue with a recombinant DNA or “transforming DNA” is hereinafter referred to as “foreign DNA” comprising one or more “transgenes”. Thus, foreign DNA may comprise both recombinant DNA as well as newly introduced, rearranged DNA of the plant.

However, the term “plant DNA” in the context of the present invention will refer to DNA of the plant that is found in the same genetic locus in the corresponding wild-type plant.

The foreign DNA can be characterized by the location and the configuration at the site of

incorporation of the recombinant DNA molecule in the plant genome. The site in the plant genome where a recombinant DNA has been inserted is also referred to as the “insertion site” or “target site”. Insertion of the recombinant DNA into the plant genome can be associated with a deletion of plant DNA, referred to as “target site deletion”. A “flanking region” or “flanking sequence” as used herein refers to a sequence of at least 20 bp, preferably at least 50 bp, and up to 5000 bp of the plant genome that is located either immediately upstream of and contiguous with or immediately downstream of and contiguous with the foreign DNA. Transformation procedures leading to random integration of the foreign DNA will result in transformants with different flanking regions, which are characteristic and unique for each transformant. When the recombinant DNA is introduced into a plant through traditional crossing, its insertion site in the plant genome, or its flanking regions will generally not be changed. An “insertion region” as used herein refers to the region corresponding to the region of at least 40 bp, preferably at least 100 bp, and up to 10000 bp, encompassed by the sequence which comprises the upstream and/or the downstream flanking region of a foreign DNA in the plant genome. Taking into consideration minor differences due to mutations within a species, an insertion region will retain, upon crossing into a plant of the same species, at least 85%, preferably 90%, more preferably 95%, and most preferably 100% sequence identity with the sequence comprising the upstream and downstream flanking regions of the foreign DNA in the plant originally obtained from transformation.

An event is defined as a (artificial) genetic locus that, as a result of genetic manipulation, carries a transgene comprising at least one copy of a gene of interest. The typical allelic states of an event are the presence or absence of the foreign DNA. An event

is characterized phenotypically by the expression of the transgene. At the genetic level, an event is part of the genetic makeup of a plant. At the molecular level, an event can be characterized by the restriction map (e.g. as determined by Southern blotting), by the upstream and/or downstream flanking sequences of the transgene, the location of

5 molecular markers and/or the molecular configuration of the transgene. Usually transformation of a plant with a transforming DNA comprising at least one gene of interest leads to a multitude of events, each of which is unique.

An elite event, as used herein, is an event that is selected from a group of events, obtained by transformation with the same transforming DNA or by back-crossing with

10 plants obtained by such transformation, based on the expression and stability of the transgene(s) and its compatibility with optimal agronomic characteristics of the plant comprising it. Thus the criteria for elite event selection are one or more, preferably two or more, advantageously all of the following:

a) That the presence of the foreign DNA does not compromise other desired

15 characteristics of the plant, such as those relating to agronomic performance or commercial value;

b) That the event is characterized by a well defined molecular configuration which is stably inherited and for which appropriate tools for identity control can be developed;

20 c) That the gene(s) of interest show(s) a correct, appropriate and stable spatial and temporal phenotypic expression, both in heterozygous (or hemizygous) and homozygous condition of the event, at a commercially acceptable level in a range of

environmental conditions in which the plants carrying the event are likely to be exposed in normal agronomic use.

It is preferred that the foreign DNA is associated with a position in the plant genome that allows easy introgression into desired commercial genetic backgrounds.

5 The status of an event as an elite event is confirmed by introgression of the elite event in different relevant genetic backgrounds and observing compliance with one, two or all of the criteria e.g. a), b) and c) above.

An “elite event” thus refers to a genetic locus comprising a foreign DNA, which answers to the above-described criteria. A plant, plant material or progeny such as seeds
10 can comprise one or more elite events in its genome.

The tools developed to identify an elite event or the plant, plant material comprising an elite event, or products that comprise plant material comprising the elite event are based on the specific genomic characteristics of the elite event, such as, a specific restriction map of the genomic region comprising the foreign DNA, molecular
15 markers or the sequence of the flanking region(s) of the foreign DNA.

Once one or both of the flanking regions of the foreign DNA have been sequenced, primers and probes can be developed which specifically recognize this (these) sequence(s) in the nucleic acid (DNA or RNA) of a sample by way of a molecular biological technique. For instance a PCR method can be developed to identify the elite event in
20 biological samples (such as samples of plants, plant material or products comprising plant material). Such a PCR is based on at least two specific “primers” one recognizing a sequence within the 5’ or 3’ flanking region of the elite event and the other recognizing a sequence within the foreign DNA. The primers preferably have a sequence of between 15

and 35 nucleotides which under optimized PCR conditions “specifically recognize” a sequence within the 5’ or 3’ flanking region of the elite event and the foreign DNA of the elite event respectively, so that a specific fragment (“integration fragment”) is amplified from a nucleic acid sample comprising the elite event. This means that only the targeted integration fragment, and no other sequence in the plant genome or foreign DNA, is amplified under optimized PCR conditions.

Preferably, the integration fragment has a length of between 50 and 500 nucleotides, most preferably of between 100 and 350 nucleotides. Preferably the specific primers have a sequence that is between 80 and 100% identical to a sequence within the 5’ or 3’ flanking region of the elite event and the foreign DNA of the elite event, respectively, provided the mismatches allow specific identification of the elite event with these primers under optimized PCR conditions. The range of allowable mismatches however, can easily be determined experimentally and are known in the art.

As the sequence of the primers and their relative location in the genome are unique for the elite event, amplification of the integration fragment will occur only in biological samples comprising (the nucleic acid of) the elite event. Preferably when performing a PCR to identify the presence of GAT-ZM1 in unknown samples, a control is included of a set of primers with which a fragment within a “housekeeping gene” of the plant species of the event can be amplified. Housekeeping genes are genes that are expressed in most cell types and which are concerned with basic metabolic activities common to all cells.

Preferably, the fragment amplified from the housekeeping gene is a fragment that is larger than the amplified integration fragment. Depending on the samples to be analyzed, other controls can be included.

Standard PCR protocols are described in the art, such as in "PCR Applications Manual" (Roche Molecular Biochemicals, 2nd Edition, 1999). The optimal conditions for the PCR, including the sequence of the specific primers, is specified in a "PCR identification protocol" for each elite event. It is however understood that a number of parameters in the PCR identification protocol may need to be adjusted to specific laboratory conditions, and may be modified slightly to obtain similar results. For instance, use of a different method for preparation of DNA may require adjustment of, for instance, the amount of primers, polymerase and annealing conditions used. Similarly, the selection of other primers may dictate other optimal conditions for the PCR identification protocol. These adjustments will however be apparent to a person skilled in the art, and are furthermore detailed in current PCR application manuals such as the one cited above.

Alternatively, specific primers can be used to amplify an integration fragment that can be used as a "specific probe" for identifying GAT-ZM1 in biological samples. Contacting nucleic acid of a biological sample, with the probe, under conditions which allow hybridization of the probe with its corresponding fragment in the nucleic acid, results in the formation of a nucleic acid/probe hybrid. The formation of this hybrid can be detected (e.g. labeling of the nucleic acid or probe), whereby the formation of this hybrid indicates the presence of GAT-ZM1. Such identification methods based on hybridization with a specific probe (either on a solid phase carrier or in solution) have been described in the art. The specific probe is preferably a sequence which, under optimized conditions, hybridizes specifically to a region within the 5' or 3' flanking region of the elite event and preferably also comprising part of the foreign DNA contiguous therewith (hereinafter referred to as "specific region"). Preferably, the specific probe

comprises a sequence of between 50 and 500 bp, preferably of 100 to 350 bp which is at least 80%, preferably between 80 and 85%, more preferably between 85 and 90%, especially preferably between 90 and 95%, most preferably between 95% and 100% identical (or complementary) to the nucleotide sequence of a specific region. Preferably, the specific probe will comprise a sequence of about 15 to about 100 contiguous nucleotides identical (or complementary) to a specific region of the elite event.

A “kit” as used herein refers to a set of reagents for the purpose of performing the method of the invention, more particularly, the identification of the elite event GAT-ZM1 in biological samples. More particularly, a preferred embodiment of the kit of the invention comprises at least one or two specific primers, as described above. Optionally, the kit can further comprise any other reagent described herein in the PCR identification protocol. Alternatively, according to another embodiment of this invention, the kit can comprise a specific probe, as described above, which specifically hybridizes with nucleic acid of biological samples to identify the presence of GAT-ZM1 therein. Optionally, the kit can further comprise any other reagent (such as but not limited to hybridizing buffer, label) for identification of GAT-ZM1 in biological samples, using the specific probe.

The kit of the invention can be used, and its components can be specifically adjusted, for purposes of quality control (e.g., purity of seed lots), detection of the elite event in plant material or material comprising or derived from plant material, such as but not limited to food or feed products.

As used herein, “sequence identity” with regard to nucleotide sequences (DNA or RNA), refers to the number of positions with identical nucleotides divided by the number of nucleotides in the shorter of the two sequences. The alignment of the two nucleotide

sequences is performed by the Wilbur and Lipmann algorithm (Wilbur and Lipmann, 1983, Proc. Natl. Acad. Sci. USA 80:726) using a window-size of 20 nucleotides, a word length of 4 nucleotides, and a gap penalty of 4. Computer-assisted analysis and interpretation of sequence data, including sequence alignment as described above, can, e.g., be conveniently performed using the programs of the Intelligenetics™ Suite (Intelligenetics Inc., CA) or the sequence analysis software package of the Genetics Computer Group (GCG, University of Wisconsin Biotechnology center). Sequences are indicated as “essentially similar” when they have a sequence identity of at least about 75%, particularly at least about 80%, more particularly at least about 85%, quite particularly about 90%, especially about 95%, more especially about 100%. When RNA sequences are said to be essentially similar or have a certain degree of sequence identity with DNA sequences, thymidine (T) in the DNA sequence is considered equal to uracil (U) in the RNA sequence. “Complementary to” as used herein refers to the complementarity between the A and T (U), and G and C nucleotides in nucleotide sequences.

The term “primer” as used herein encompasses any nucleic acid that is capable of priming the synthesis of a nascent nucleic acid in a template-dependent process, such as PCR. Typically, primers are oligonucleotides from 10 to 30 basepairs, but longer sequences can be employed. Primers may be provided in double-stranded form, though the single-stranded form is preferred. Probes can be used as primers, but are designed to bind to the target DNA or RNA and need not be used in an amplification process.

The term “recognizing” as used herein when referring to specific primers, refers to the fact that the specific primers specifically hybridize to a nucleic acid sequence in the

elite event under the conditions set forth in the method (such as the conditions of the PCR identification protocol), whereby the specificity is determined by the presence of positive and negative controls.

The term “hybridizing” as used herein when referring to specific probes, refer to the fact that the probe binds to a specific region in the nucleic acid sequence of the elite event under standard stringency conditions. Standard stringency conditions as used herein refers to the conditions for hybridization described herein or to the conventional hybridizing conditions as described by Sambrook et al. (1989) (Molecular Cloning: A Laboratory Manual, Second Edition, Cold Spring Harbor Laboratory Press, NY) which for instance can comprise the following steps: 1) immobilizing plant genomic DNA fragments on a filter, 2) prehybridizing the filter for 1 to 2 hours at 42°C in 50% formamide, 5 X SSPE, 2 X Denhardt’s reagent and 0.1% SDS, or for 1 to 2 hours at 68°C in 6 X SSC, 2 X Denhardt’s reagent and 0.1% SDS, 3) adding the hybridization probe which has been labeled, 4) incubating for 16 to 24 hours, 5) washing the filter for 20 min. at room temperature in 1X SSC, 0.1 %SDS, 6) washing the filter three times for 20 min. each at 68°C in 0.2 X SSC, 0.1 %SDS, and 7) exposing the filter for 24 to 48 hours to X-ray film at -70°C with an intensifying screen.

As used in herein, a biological samples is a sample of a plant, plant material or products comprising plant material. The term “plant” is intended to encompass corn (*Zea mays*) plant tissues, at any stage of maturity, as well as any cells, tissues, or organs taken from or derived from any such plant, including without limitation, any seeds, leaves, stems, flowers, roots, single cells, gametes, cell cultures, tissue cultures or protoplasts. “Plant material”, as used herein refers to material which is obtained or derived from a

plant. Products comprising plant material relate to food, feed or other products that are produced using plant material or can be contaminated by plant material. It is understood that, in the context of the present invention, such biological samples are tested for the presence of nucleic acids specific for GAT-ZM1, implying the presence of nucleic acids in the samples. Thus the methods referred to herein for identifying elite event GAT-ZM1 in biological samples, relate to the identification in biological samples of nucleic acids which comprise the elite event.

As used herein “comprising” is to be interpreted as specifying the presence of the stated features, integers, steps, reagents or components as referred to, but does not preclude the presence or addition of one or more features, integers, steps or components, or groups thereof. Thus, e.g., a nucleic acid or protein comprising a sequence of nucleotides or amino acids, may comprise more nucleotides or amino acids than the actually cited ones, i.e., be embedded in a larger nucleic acid or protein. A chimeric gene comprising a DNA sequence that is functionally or structurally defined, may comprise additional DNA sequences, etc.

The following examples describe the identification of the development of tools for the identification of elite event GAT-ZM1 in biological samples.

Unless otherwise stated, all recombinant DNA techniques are carried out according to standard protocols as described in Sambrook et al. (1989) *Molecular Cloning: A Laboratory Manual*, Second Edition, Cold Spring Harbor Laboratory Press, NY and in Volumes 1 and 2 of Ausubel et al. (1994) *Current Protocols in Molecular Biology*, Current Protocols, USA. Standard materials and methods for plant molecular work are

described in Plant Molecular Biology Labfax (1993) by R.D.D. Croy published by BIOS Scientific Publications Ltd (UK) and Blackwell Scientific Publications, UK.

In the description and examples, reference is made to the following sequences:

- SEQ ID NO: 1: sequence of the genetic elements of vector pUC/Ac
- 5 SEQ ID NO: 2: primer MDB286
- SEQ ID NO: 3: primer MDB391
- SEQ ID NO: 4: primer MDB411
- SEQ ID NO: 5: primer MDB420
- SEQ ID NO: 6: nucleotide sequence comprising a 5' flanking region of GAT-ZM1
- 10 SEQ ID NO: 7: primer MDB439
- SEQ ID NO: 8: primer VDS44
- SEQ ID NO: 9: primer MDB522
- SEQ ID NO: 10: nucleotide sequence comprising a 3' flanking region of GAT-ZM1
- SEQ ID NO: 11: primer COR17
- 15 SEQ ID NO: 12: primer COR18
- SEQ ID NO: 13: primer COR15
- SEQ ID NO: 14: primer COR16

Example 1

Identification of the flanking regions of elite event GAT-ZM1

20 Herbicide-resistant corn was developed by transformation of corn with the pUC/Ac vector comprising the coding sequence of a *pat* gene encoding the enzyme phosphinothricin-acetyl-transferase, under the control of the constitutive 35S promoter from Cauliflower Mosaic virus (CaMV). A detailed description of the genetic elements of

pUC/Ac is provided in Table 1. The nucleotide sequence of the genetic elements of pUC/Ac is provided in SEQ ID NO: 1.

Table: genetic elements of vector pUC/Ac

Nuc. No:	Genetic element
412 – 618	35S terminator from CaMV from the vector pDH51. Pietrzak et al. (1986) Nucl. Acids Res. 14:5857-5868
619 – 636	Synthetic polylinker sequences
637 – 1188	Synthetic pat gene (amino acid sequence from <i>Streptomyces viridochromogenes</i>). European Patent 275957 B1
1189 – 1216	Synthetic polylinker sequences
1217 – 1746	35S promoter from CaMV from the vector pDH51. Pietrzak et al. (1986).
1747 – 411	Sequence of the vector pUC18, including the β -lactamase gene (pos. 2923 – 3783) and the origin of replication at pos. 2164. Yanisch-Perron et al. (1985) Gene 33:103-119

Elite event GAT-ZM1 was selected based on an extensive selection procedure based on good expression and stability of the herbicide resistance gene and its compatibility with optimal agronomic characteristics.

The sequence of the regions flanking the foreign DNA in the GAT-ZM1 event was determined using the thermal asymmetric interlaced (TAIL-) PCR method described by Liu et al. ((1995) Plant J. 8:457-463). This method utilizes three nested primers in successive reactions together with a shorter arbitrary degenerate primer so that the relative amplification efficiencies of specific and non-specific products can be thermally controlled. The specific primers were selected for annealing to the border of the foreign DNA and based on their annealing conditions. A small amount (5 μ l) of unpurified secondary and tertiary PCR products were analyzed on a 1% agarose gel. The tertiary PCR product was used for preparative amplification, purified and sequenced on an automated sequencer using the DyeDeoxy Terminator cycle kit.

1.1. Right (5') flanking region

The primers used were:

	Sequence (5' → 3')	Position in pUC/Ac
Degenerate primer MDB286	NgT.CgA.SWg.ANA.WgA.A (SEQ ID NO: 2)	-
Primary TAIL MDB391	Tgg.ATA.CAA.gCA.Tgg.Tgg.ATg.g (SEQ ID NO: 3)	715←736
Secondary TAIL MDB411	Agg.CAT.gCC.gCT.gAA.ATC.ACC (SEQ ID NO: 4)	606←626
Tertiary TAIL MDB420	GgT.TTC.gCT.CAT.gTg.TTg.AgC (SEQ ID NO: 5)	507←527

Whereby: N = A,C,T or g; S = C or g; W = A or T

The fragment amplified using MDB286-MDB420 was ca. 1100 bp, the complete sequence of that was determined (SEQ ID NO: 6). The sequence between nucleotide 1 and 341 corresponds to plant DNA, while the sequence between nucleotide 342 and 1041 corresponds to T-DNA.

1.2. Left (3') flanking region

The primers used were:

	Sequence (5' → 3')	Position in pUC/Ac
Degenerate primer MDB286	NgT.CgA.SWg.ANA.WgA.A (SEQ ID NO: 2)	-
Primary TAIL MDB439	CTC.ATg.gTT.ATg.gCA.gCA.CTg.C (SEQUENCE ID NO: 7)	3401-3422
Secondary TAIL VDS44	CTg.TCA.TgC.CAT.CCG.TAA.gAT.gC (SEQUENCE ID NO: 8)	3435-3457
Tertiary TAIL MDB522	TgC.TTT.gAA.gAC.gTg.gTT.gg (SEQUENCE ID NO: 9)	1326-1345

Whereby: N = A,C,T or g; S = C or g; W = A or T

The fragment amplified using MDB286-MDB522 was ca. 450 bp, the complete sequence of which was determined (SEQ ID NO: 10). The sequence between nucleotide 1

and 342 corresponds to T-DNA, while the sequence between nucleotide 343 and 484 corresponds to plant DNA.

Example 2 Development of a Polymerase Chain reaction identification protocol

2.1. Primers

Specific primers were developed which recognize sequences within the elite event. More particularly, a primer was developed which recognizes a sequence within the 5' flanking region of GAT-ZM1. A second primer was then selected within the sequence of the foreign DNA so that the primers span a sequence of about 200 bp. The following primers were found to give particularly clear and reproducible results in a PCR reaction on GAT-ZM1 DNA:

Name	Sequence	SEQ ID NO:	Target:
COR17:	5'-ggg.TgA.gCT.CgA.ATg.TTg.TTC.T-3'	11	plant DNA
COR18:	5'-TCT.TAg.ACg.TCA.ggT.ggC.ACT.T-3'	12	T-DNA

Primers targeting an endogenous sequence are preferably included in the PCR cocktail. These primers serve as an internal control in unknown samples and in the DNA positive control. A positive result with the endogenous primer-pair demonstrates that there is ample DNA of adequate quality in the genomic DNA preparation for a PCR product to be generated. The endogenous primers were selected to recognize a housekeeping gene in *Zea mays*:

Name	Sequence	SEQ ID NO	Located in
COR15:	5'-AgC.gTC.AAg.gAT.CAT.Tgg.TgT.C-3'	13	<i>Zea Mays</i> alcohol dehydrogenase 1 gene (X04050)
COR16:	5'-ggC.CAA.gTT.CAg.CAT.AAg.CTg.T-3'	14	<i>Zea Mays</i> alcohol dehydrogenase 1 gene (X04050)

2.2. Amplified fragments

The expected amplified fragments in the PCR reaction are:

For primer pair COR15-COR16: 513bp (endogenous control)

For primer pair COR17-COR18: 202bp (GAT-ZM1 elite Event)

2.3. Template DNA

Template DNA was prepared from a leaf punch according to Edwards et al. ((1991) Nucl. Acids Res. 19:1349). When using DNA prepared with other methods, a test run utilizing different amounts of template should be done. Usually 50 ng of genomic template DNA yields the best results.

2.4. Assigned positive and negative controls

To avoid false positives or negatives, it was determined that the following positive and negative controls should be included in a PCR run:

- Master Mix control (DNA negative control). This is a PCR in which no DNA is added to the reaction. When the expected result, no PCR products, is observed this indicates that the PCR cocktail was not contaminated with target DNA.
- A DNA positive control (genomic DNA sample known to contain the transgenic sequences). Successful amplification of this positive control demonstrates that the PCR was run under conditions that allow for the amplification of target sequences.
- A wildtype DNA control. This is a PCR in which the template DNA provided is genomic DNA prepared from a non-transgenic plant. When the expected result, no amplification of a transgene PCR product but amplification of the

endogenous PCR product, is observed this indicates that there is no detectable transgene background amplification in a genomic DNA sample.

2.5. PCR conditions

Optimal results were obtained under the following conditions:

- the PCR mix for 25µl reactions contains:

- 2.5 µl template DNA
 - 2.5 µl 10x Amplification Buffer (supplied with Taq polymerase)
 - 0.5 µl 10 mM dNTP's
 - 0.5 µl COR17 (10pmoles/µl)
 - 0.5 µl COR18 (10pmoles/µl)
 - 0.25 µl COR15 (10pmoles/µl)
 - 0.25 µl COR16 (10pmoles/µl)
 - 0.1 µl Taq DNA polymerase (5 units/µl)
 - water up to 25 µl

- the thermocycling profile to be followed for optimal results is the following:

- 4 min. at 95°C

- Followed by:

- 1 min. at 95°C

- 1 min. at 57°C

- 2 min. at 72°C

- For 5 cycles

- Followed by:

- 30 sec. at 92°C

- 30 sec. at 57°C

- 1 min. at 72°C

- For 25 cycles

Followed by:

5 minutes at 72°C

2.6. Agarose gel analysis

To optimally visualize the results of the PCR it was determined that between 10
5 and 20µl of the PCR samples should be applied on a 1.5% agarose gel (Tris-borate buffer)
with an appropriate molecular weight marker (e.g. 100bp ladder PHARMACIA).

2.7. Validation of the results

It was determined that data from transgenic plant DNA samples within a single
PCR run and a single PCR cocktail should not be acceptable unless 1) the DNA positive
10 control shows the expected PCR products (transgenic and endogenous fragments), 2) the
DNA negative control is negative for PCR amplification (no fragments) and 3) the wild-
type DNA control shows the expected result (endogenous fragment amplification).

When following the PCR Identification Protocol for GAT-ZM1 as described
above, lanes showing visible amounts of the transgenic and endogenous PCR products of
the expected sizes, indicate that the corresponding plant from which the genomic template
15 DNA was prepared, has inherited the GAT-ZM1 elite event. Lanes not showing visible
amounts of either of the transgenic PCR products and showing visible amounts of the
endogenous PCR product, indicate that the corresponding plant from which the genomic
template DNA was prepared, does not comprise the elite event. Lanes not showing visible
20 amounts of the endogenous and transgenic PCR products, indicate that the quality and/or
quantity of the genomic DNA didn't allow for a PCR product to be generated. These
plants cannot be scored. The genomic DNA preparation should be repeated and a new
PCR run, with the appropriate controls, has to be performed.

2.8. Use of discriminating PCR protocol to identify GAT-ZM1

Before attempting to screen unknowns, a test run, with all appropriate controls, has to be performed. The developed protocol might require optimization for components that may differ between labs (template DNA preparation, Taq DNA polymerase, quality of the primers, dNTP's, thermocycler, etc.).

Amplification of the endogenous sequence plays a key role in the protocol. One has to attain PCR and thermocycling conditions that amplify equimolar quantities of both the endogenous and transgenic sequence in a known transgenic genomic DNA template. Whenever the targeted endogenous fragment is not amplified or whenever the targeted sequences are not amplified with the same ethidium bromide staining intensities, as judged by agarose gel electrophoresis, optimization of the PCR conditions may be required.

Zea mays leaf material from a number of plants, some of which comprising GAT-ZM1 were tested according to the above-described protocol. Samples from elite event GAT-ZM1 and from *Zea mays* wild-type were taken as positive and negative controls, respectively.

Figure 1 illustrates the result obtained with the elite event PCR identification protocol for GAT-ZM1 on a number of corn plant samples (lanes 1 to 14). The samples in lanes 1, 2, 5, 6, 7, 8, 11, 13, and 14 were found to contain the elite event as the 202 bp band is detected, while the samples in lanes 4, 9, 10, and 12 do not comprise GAT-ZM1. Lane 3 indicates a PCR failure, as the control band is not detected. Lanes 19 and 20 represent GAT-ZM1 positive control samples, lanes 20 and 22 represent non-transgenic *Zea mays* controls; lane 23 represents the negative control (water) sample, and lane 24 the Molecular Weight Marker (100 bp).

Example 3
Use of a specific integration fragment as a probe for
detection of material comprising GAT-ZM1

A specific integration fragment of GAT-ZM1 is obtained by PCR amplification
5 using specific primers COR17 (SEQ ID NO: 11) and COR18 (SEQ ID NO: 12) or by
chemical synthesis and is labeled. This integration fragment is used as a specific probe for
the detection of GAT-ZM1 in biological samples. Nucleic acid is extracted from the
samples according to standard procedures. This nucleic acid is then contacted with the
specific probe under hybridization conditions that are optimized to allow formation of a
10 hybrid. The formation of the hybrid is then detected to indicate the presence of GAT-ZM1
nucleic acid in the sample. Optionally, the nucleic acid in the samples is amplified using
the specific primers prior to contact with the specific probe. Alternatively, the nucleic
acid is labeled prior to contact with the specific probe instead of the integration fragment.
Optionally, the specific probe is attached to a solid carrier (such as, but not limited to a
15 filter, strip or beads), prior to contact with the samples.